

		Minitom
Capacity	Height x Length	30 x 40 mm (1.2" x 1.6")
	Diameter	40 mm (1.6")
	Maximum cutting force	0.34 N
Cut-off wheel	Diameter	100-127 mm (4"-5")
	Arbor diameter	12.7 mm (0.5")
Cut-off wheel motor	Rotational speed	100-420 rpm
Recirculation cooling tank	Volume	250 ml (0.07 gal)
Software and electronics	Controls	Touch pad
	Display	N/A
Safety standards	See the Declaration of Conformity	
REACH	For information about REACH, contact your local Struers office.	
Operating environment	Surrounding temperature	5-40 °C (41 - 131 °F)
	Humidity	< 95 % RH non-condensing
	Noise level	L(Pa)=58 dBa(A)
Power supply	Voltage/frequency	100-250 AVC 50-60 Hz
	Power inlet	1L + (N) + PE
	Power S1	N/A
	Current, nominal load	3.2 A
	Current, max.	4.8 A
	Ampere rating of the largest motor or load	N/A
Dimensions and weight	Width	27.5 cm (10.8")
	Depth (main body)	28 cm (11")
	Depth (with handle)	41 cm (16")
	Height (main body)	20.5 cm (8.1")
	Height (with handle)	25.6 cm (10.1")
	Weight	8 kg (17.6 lbs)